

# A New Generation of Rapid Assays: Detection of HIV and Syphilis Simultaneously

Niel T. Constantine, Ph.D., MT (ASCP)

Professor of Pathology

University of Maryland School of Medicine

Institute of Human Virology

Director, Laboratory of Viral Diagnostics

March 27, 2019

No conflict of Interest




# Rapid Tests: Advantages

- Rapid results
- Appropriate for Resource-limited Countries
  - No stable electricity required
  - No cold chain required for shipment
  - Robust – storage up to 30C
- Simplified training of staff & “no training necessary”
- Only a single finger-stick or oral fluid required

# History


# Multiplex Rapid Tests

# *Platform Device Other Applications*




**Universal HealthWatch  
Diagnostics Corp.** 

# Rapid Confirmatory Test



# Determine® HIV-1/2 Ag/Ab Combo



# Multiplex Rapid Assays

- DPP HIV-Syphilis System (Chembio Diagnostic Systems)
- INSTI Multiplex HIV-1/2 Syphilis Ab Test
- SD Bioline HIV/Syphilis Duo (Standard Diagnostics, Alere)
- First Response HIV1+2/Syphilis COMBO Card Test (Premier Medical)
- First Response HIV1+2/HCV COMBO Card Test (Premier Medical)
- SD Bioline Roto/Adeno (Standard Diagnostics)
- HIV/HCV/HBsAg/Syphilis Combo Rapid Test (Healgen)
- Others:



# World Health Organization

## Recommendations for Use of Dual HIV/Syphilis Rapid Diagnostic Tests (2017)

- First Dual HIV/Syphilis test listed on WHO list of Prequalified *in vitro* diagnostic products in 2015
- Final sensitivity/specificity for HIV: 100 / 99.5
- Final sensitivity/specificity for Syphilis: 87 / 99.5
- WHO recommends HIV and syphilis testing for woman at least once, preferably in the first trimester, and WHO states that “the use of the Dual HIV/Syphilis RTD would facilitate implementation of this recommendation”.

Clinical Trial

Evaluation of the  
DPP HIV-Syphilis System

# DPP® HIV-Syphilis System


## 2017-2018




**Manufacturer: Chembio Diagnostics Systems, Inc.**

# The DPP HIV-Syphilis System


10-30 minutes total time



# Negative for HIV and Syphilis



# Positive for HIV and Syphilis



# DPP Micro reader



# DPP® Micro Reader



DPP Micro Reader



Test Device Holder



DPP Micro Reader With Test  
Device Holder and Test Device

# Procedures

- Gather subject information
  - Hard copy (eligibility, contact information)
  - Directly enter demographics and history into interactive database
- Collect 4 tubes of blood
  - Centrifuge within 2 hours
- Collect fingerstick specimen and test on DPP
- Test whole blood and plasma on DPP
- Ship tubes of blood to reference laboratory

# Samples Tested on Site

---

- Whole Blood
- Fingerstick
- Plasma

# Reference Testing at Central Laboratory

- HIV Ag/Ab 4<sup>th</sup> generation test (Abbott)
  - Multispot Rapid HIV-1/2 differentiation test or Geenius HIV-1/2 (Bio-Rad)
  - HIV RNA if two tests were discordant (Gen-Probe Aptima)
- 
- Syphilis Screen: Bio-Rad Bioplex 2200 Syphilis IgG EIA (treponemal)
  - Syphilis RPR (non-treponemal) with titer
  - Syphilis TP PA test (treponemal)

# Population

- 2 sites:
  - University of Maryland, Baltimore, Maryland
  - Fanno Creek Clinic, Portland, Oregon

- Total Sample Size with all matrices: 803

UMB: 324

Fanno Creek Clinic: 479

HIV Known positives

Syphilis Known Positives

Known Co-infected

High risk unknown status

Low risk

# Results

## Reference Tests

HIV Positive:	351
Syphilis Positive:	119
Co-infected:	91
Negative:	
HIV	427
Syphilis	672

# Results HIV

## DPP HIV-Syphilis System: Summary of Overall HIV Sensitivity and Specificity (n=778)

Sample Matrix	HIV Sensitivity (95% CI)	HIV Specificity (95% CI)
Fingerstick Blood	$348/350^1 = 99.4\% (97.9\% - 99.8\%)$	$421/424^2 = 99.3\% (97.9\% - 99.8\%)$
Venous Whole Blood	$350/351 = 99.7\% (98.4\% - 99.9\%)$	$425/427 = 99.5\% (98.3\% - 99.9\%)$
Plasma	$350/351 = 99.7\% (98.4\% - 99.9\%)$	$424/427 = 99.3\% (98.0\% - 99.8\%)$

<sup>1</sup>Fingerstick Whole Blood Results not collected from one (1) subject.

<sup>2</sup>Fingerstick Whole Blood Results not collected from three (3) subjects.

# Results Syphilis

## DPP HIV-Syphilis System: Summary of Overall *T. pallidum* Sensitivity and Specificity (n=791)

Sample Matrix	Syphilis Sensitivity (95% CI)	Syphilis Specificity (95% CI)
Fingerstick Blood	$111/119 = 93.3\% (87.3\% - 96.6\%)$	$646/668^1 = 96.7\% (95.1\% - 97.8\%)$
Venous Whole Blood	$115/119 = 96.6\% (91.7\% - 98.7\%)$	$646/672 = 96.1\% (94.4\% - 97.3\%)$
Plasma	$116/119 = 97.5\% (92.9\% - 99.1\%)$	$647/672 = 96.3\% (94.6\% - 97.5\%)$

<sup>1</sup>Fingerstick Whole Blood Results not collected from four (4) subjects.

# False Negatives: Seroreversion?

## Comparative Percent T. pallidum Reactivity in Plasma with Previously Diagnosed Syphilis Subjects

DPP® HIV-Syphilis System  
(Rapid Test)


85/107 = 79.4%  
95% CI: 70.8–86.0%

Bio-Rad Bioplex 2200  
(Automated Multi-test Platform)

85/107 = 79.4%  
95% CI:70.8–86.0%

# SYPHILIS REACTIVE SAMPLES

## READER VALUE VS RPR TITER



# Future Applications of the DPP Reader

