



Department
of Health

Wadsworth
Center

Evaluation of the Aptima HCV Quant Dx Assay for Qualitative and Quantitative RNA Detection Using Dried Blood Spots

Renee Hallack

Wadsworth Center

New York State Department of Health

renee.hallack@health.ny.gov

(518)474-4882



Hepatitis C virus (HCV) Molecular Testing in the Bloodborne Virus Lab (BVL)

15 Years

- HCV Genotyping with Phylogenetic Analysis
 - Outbreak Investigations
 - Surveillance

10 Years

- Confirmatory HCV RNA Testing
- HCV Viral Load Testing

3 Years

- Dried blood spots (DBS) for Confirmatory HCV RNA Testing

This Year

- Seeking Approval to use DBS for HCV Viral Load Testing



- ✓ Random Access
- ✓ High Throughput
- ✓ Test Results in 3 Hours

Hologic's Aptima HCV Quant DX Assay is FDA Approved

- Real-Time Transcription Mediated Amplification Test
- Fully Automated Panther System
- Limit of Detection (LOD) of 3.9 IU/mL
- Lower Limit of Quantitation (LLOQ) of 10 IU/mL
- Fresh or Frozen Plasma or Serum Specimens




Department
of Health



Recommended Testing Sequence for Identifying Current Hepatitis C Virus (HCV) Infection





U.S. Department of Health and Human Services
Centers for Disease Control and Prevention



Department of Health





They are simply collected by fingerstick onto filter paper.

**DBS:
An
Alternative
to Plasma**



No phlebotomist, centrifuges or cold chain needed.

*Photos from the New York State Department of Health's HCV Rapid Testing Program's Training for HCV RNA Testing via Dried Blood Spot



Department of Health



All DBS for validation were made in the BVL Not from fingerstick

**DBS
Used for
Validation**

2 Types:

- From Venipuncture
HCV positive EDTA-whole blood
Spotted prior to separating plasma

- Contrived
HCV negative whole blood
OR red blood cells spiked
with HCV positive plasma



Department
of Health



Clinical Laboratory Evaluation Program (CLEP) Modified Specimen Type for FDA Approved Assay Guidelines- Validation MUST Include:

1. Limit of Detection (LOD)
2. Lower Limit of Quantitation (LLOQ)
3. Upper Limit of Quantitation (ULOQ)
4. Stability at room temperature and 37°C, within 1 Log of LLOQ
5. Intra-Assay Reproducibility
6. Inter-Assay Reproducibility
7. Accuracy

Note: For 2-7 All viral load results must be compared to plasma

How can we compare 2 6-mm punches of a DBS
to 700 uL of plasma?



Department
of Health



Determining DBS Conversion Factor



Journal of Clinical Virology

Volume 112, March 2019, Pages 40-44



Performance evaluation of the Hologic Aptima HCV Quant Dx assay for detection of HCV RNA from dried blood spots

Beth Catlett^{a, b, c, d}, Alex Carrera^a, Mitchell Starr^a, Tanya L Applegate^b, Peter Lowe^c, Jason Grebely^{b, 1}, H. Philip Cunningham^{a, 1}

- $\text{volume of ATM } (\mu\text{l}) / [\text{DBS volume } (\mu\text{l}), \times (1 - .45)]$
 - $(1000/70 \times .55 = 25.97)$



Department
of Health



Volume of Buffer/ DBS Volume x 1-Hematocrit

Volume of (DBS Elution Buffer)= 750 uL

DBS Volume=2 x 10 uL Punches= 20 uL

1-Hemarocrit= 1-0.45= 0.55

**Determining
OUR DBS
Conversion
Factor**

750 (μl) / [20 (μl), x 0.55]

750 (μl) / 11= 68.2



Department
of Health



Limit of Detection (LOD) Goals:

- ✓ Duplicate testing on multiple runs
- ✓ $\geq 95\%$ detected (qualitative)

DBS LOD Shown as percent detected for two 6-mm Punches

IU/mL	# Replicates tested	# Detected	% Detected
1000	7	7	100%
500	7	7	100%
250	7	7	100%
125	7	1	14%
62.5	7	3	43%
31.25	7	0	0%

Result: We can detect down to
2.4 Log₁₀ IU/mL (250 IU/mL)
This is our LOD



Department
of Health



Lower Limit of Quantitation (LLOQ) Goals:

- ✓ 20 replicates on multiple runs
- ✓ $\geq 95\%$ quantitated within 1 Log of plasma

Dilution 

Log ₁₀ IU/mL	3.12	3.14	3.26	3.34	3.42	3.57	3.61
# tested	20	20	20	20	20	20	20
# quantitated	19	20	20	19	20	20	20
% quantitated	95	100	100	95	100	100	100
% with VL diff ≤ 1	0	0	0	0	0	0	0
Log difference range	-0.35-0.13	-0.38-0.03	-0.25-0.42	-0.31-0.39	-0.10-0.44	-0.15-0.42	-0.21-0.41

Result: We can quantitate down to

3.12 Log₁₀ IU/mL (1304 IU/mL)

This is our LLOQ




Department
of Health



Upper Limit of Quantitation (ULOQ) Goals:

- ✓ Quantitated result within 1 log of plasma
- ✓ 20 replicates on multiple runs
- ✓ Low SD between samples

Viral Load Results for 2 Plasma (Mean) and 20 DBS
Replicates Tested to Determine ULOQ



Run	N	[Mean] IU/mL	[Mean] log ₁₀ IU/mL	Intra Run SD
1	10	37,033,041	7.61	0.06
2	7	44,855,939	7.65	0.07
3	3	49,457,526	7.69	0.05
Total	20		Inter Run SD (Total SD)	0.08

7.65 Log₁₀ IU/ml

The highest viral load we can accurately quantitate



Department of Health



Stability Question: At what condition does the quantitative value become unreliable?

Plasma VL 3.61	Storage Condition	
Sample	Ambient	37°C
3dayRep 1	3.41	3.57
3dayRep 2	3.49	3.44
3dayRep 3	3.4	3.62
6dayRep 1	3.49	3.38
6dayRep 2	3.51	3.25
6dayRep 3	3.52	3.47
8dayRep 1	3.51	3.45
8dayRep 2	3.71	3.4
8dayRep 3	3.27	3.4
14dayRep 1	3.44	3.66
14dayRep 2	Not done	3.53
14dayRep 3	3.33	3.38
16dayRep 1	3.54	3.31
16dayRep 2	3.62	3.52
16dayRep 3	3.75	3.75
Mean	3.5	3.47
SD	0.13	0.14
% CV	3.81	3.9

We still don't know the answer

We do know that:
Little variability exists between our conditions across all timepoints and when compared to plasma



Reproducibility:

Intra-Assay: 3 replicates of 5 specimens

Inter-Assay: 3 runs of 4 triplicate specimens

Sample	Dilution 1	Dilution 2	Dilution 3	Dilution 4	Dilution 5
Intra-Assay Reproducibility Run 1					
Mean	3.35	3.83	4.27	4.48	4.71
SD	0.12	0.09	0.07	0.09	0.04
% CV	3.59	2.27	1.54	2.07	0.90
Intra-Assay Reproducibility Run 2					
Mean	3.28	3.71	4.12	4.49	4.85
SD	0.10	0.10	0.04	0.04	0.10
% CV	3.07	2.71	0.95	0.84	2.11
Intra-Assay Reproducibility Run 3					
Mean	Sample Exhausted	3.78	4.11	4.32	4.72
SD		0.10	0.05	0.06	0.06
% CV		2.63	1.22	1.44	1.18
Inter-Assay Reproducibility Run 1- Run 3					
Mean	Not Done	3.77	4.17	4.43	4.76
SD		0.06	0.09	0.10	0.08
% CV		1.55	2.07	2.17	1.58

[SD Range] 0.04-0.12
[%CV] 0.84- 3.59

Little variability within runs and between runs




Blinded Accuracy

One plasma sample
viral load (VL)=15
(<DBS LOD)

	Quantitated Plasma	Negative Plasma	Total
Quantitated DBS	29	0	29
Negative DBS	1	11	12
Total	30	11	41

R^2 is acceptable

The lowest quantitated
DBS VL= 748



There is a strong correlation between plasma and DBS




Department
of Health



Linear Regression

DBS Adjusted VL vs Plasma VL
N=78 $R^2 = 0.9449$




Predicted strong,
positive correlation
between plasma
and DBS holds true



Department
of Health



Bland Altman Analysis: 78 Paired Specimens



Ideally, bias should be close to 0 if producing similar results: **-0.03**
≥95% of points should be between the Limits of Agreement (LOA): **96.9%**

There is a high degree of agreement between plasma and DBS

The assay is accurate and highly reproducible when using the DBS conversion factor of 68.2 for DBS:

- ✓ stored up to 15 days between room temperature and 37°C
- ✓ viral loads between 3.12 and 7.98 Log₁₀ IU/mL

**Should we
use DBS for
viral load
testing?**

**Yes
If
venipuncture
is not an
option**



**Department
of Health**



Acknowledgements

- The New York State Department of Health's AIDS Institute
- The AIDS Institute's HCV Rapid Testing Program
- New York State Rapid Testing Sites
- The Wadsworth Center's Bloodborne Virus Laboratory



Department
of Health

